

Abundant Clinical Applications




Placenta B Mode



Umbilical Cord



Fetal Face 4D Mode




Panoramic Imaging



IMT Auto Measurement



Liver Contrast-enhanced



Hepatic Vein C Mode



First Hepatic Portal



Fetal Leg



Wuhan Zoncare Bio-medical Electronics Co.,Ltd.

Address: Zoncare Building, No. 380, High-tech 2nd Road, Eastlake High-tech Development Zone, Wuhan, Hubei 430206 P. R. China.
Tel: +86-27-8777 0203 Fax: +86-27-8777 0581

Email: info@zoncare.com <http://www.zoncare.com>

© 2020 Wuhan Zoncare Bio-Medical Electronics Co.,Ltd. All rights reserved. Specifications subject to changes without prior notice.



ZONCARE-N7

Full Digital Color Doppler Ultrasonic Diagnostic System



<http://www.zoncare.com>



Features

- 19-inch high resolution medical LED monitor
- 10.4-inch touch screen
- Four activated transducer connectors
- Standard configured gel heating cup
- Carotid IMT auto measurement
- Optional 3D/4D imaging
- Dual/Triplex synchronous technology
- Spectrum full screen display
- Support Multi-parameter presets
- DICOM 3.0 supported
- Support multiple data transmission
- Support HDMI output interface

Advanced Imaging Technology



Multiple Imaging Modes

- Color Flow Imaging (CFM)
- Pulse Wave Doppler Imaging(PW)
- Power Doppler Imaging(PDI)
- Directional Power Doppler Imaging(DPDI)
- Continuous Wave Doppler(CW)
- High Pulse Repetition Frequency(HPRF)
- Tissue Doppler Imaging(TDI)
- Color Motion Mode(CMM)
- Anatomic Motion Mode(AMM)
- Panoramic Imaging(Wfov)
- Elastic Imaging(E)
- 3D/4D mode

